

Sterile Respirator Solutions

Protection for your worker. Assurance for your product.

Bullard sterile respirator solutions are designed to help you maintain your sterile environments. You'll find assurance in the pre-assembled, gamma irradiated, double-bagged selection of sterile hoods and facepieces. Pair Bullard sterile hoods and facepieces with a Bullard Powered Air-Purifying Respirator (PAPR) for a complete package of protection, comfort and value for your workers.



Superior Worker Protection

- Full hoods offer 1,000 Assigned Protection Factor (APF)
- Loose-fitting facepieces provide 150% more protection than N95 masks
- Patented inflatable neck cuff helps ensure positive pressure
- DuPont™ Tychem® QC construction offers superior protection

For use with the industry-leading Bullard Powered Air Purifying Respirators (PAPRs)

- Intrinsically safe (IS) protection available
- Audible low voltage alarms help protect workers
- 33% more airflow delivered than NIOSH minimum
- Up to ten hours of run time provides strong performance throughout a shift

Sterile Product Confidence

- Gamma irradiated
- Pre-assembled
- Double bagged
- 10^{-6} Sterility Assurance Level
- FDA-validated



Sterilization...what is it?

Sterilization is a process eliminating viable microorganisms from an item or area. The sterility of an item is measured as a probability, which is called the Sterility Assurance Level. Although sterilization significantly reduces the probability of the presence of microorganisms, the existence of microorganisms can never be reduced to zero with complete certainty.



*PAPRs mentioned in this document are not gamma irradiated, and therefore no claim is made regarding sterility assurance for this component.



20TICRSA Sterile Respirator Solution
Sterile assembly includes:

- 20TICR Double bib hood with Ratchet Suspension
- PA1BT Breathing Tube
- Fully Assembled
- Packaged Five per box.
- Minimum Quantity: 100 assemblies (20 boxes)


RT3SA Sterile Respirator Solution
Sterile assembly includes:

- RT3 Double bib hood
- PA1BT Breathing tube
- Fully assembled
- Packaged five per box
- Minimum quantity: 100 assemblies (20 boxes)


20LF2L & 20LF2M Sterile Respirator Solution
Sterile assembly includes:

- 20LF2L (large) or 20LF2M (medium)
- PA20LFBT Breathing Tube
- Fully Assembled
- Packaged Five per box.
- Minimum Quantity: 100 assemblies (20 boxes)


20TICHSA Sterile Respirator Solution
Sterile assembly includes:

- 20TICH Double bib hood for use with hard hats
- PA1BT Breathing Tube
- Fully Assembled
- Packaged Five per box.
- Minimum Quantity: 100 assemblies (20 boxes)


Bullard

1898 Safety Way
Cynthiana, KY 41031-9303
Toll free: 877-BULLARD (285-5273)
Phone: 859-234-6616
Fax: 859-234-8987
www.bullard.com

Bullard GmbH

Lilienthalstrasse 12
53424 Remagen
Germany
Phone: +49-2642 999980
Fax: +49-2642 9999829
www.bullardextrem.com

Bullard Asia Pacific PTE. LTD.

LHK Building
701, Sims Drive, #04-03
Singapore 387383
Phone: (65) 6745-0556
Fax: (65) 6745-5176
www.bullard.com

www.bullard.com



©2009 Bullard. All rights reserved.
Tychem is a registered trademark of E.I. du Pont de Nemours and Company.